

AQ 1989

Chemical Properties


CAS No. : 61-40-5

Formula: C₁₃H₁₆ClNOS

Molecular Weight: 269.79

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	AQ 1989 is a platelet aggregation inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7066 mL	18.5329 mL	37.0659 mL
5 mM	0.7413 mL	3.7066 mL	7.4132 mL
10 mM	0.3707 mL	1.8533 mL	3.7066 mL
50 mM	0.0741 mL	0.3707 mL	0.7413 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481