


AI 3-22016

## Chemical Properties

CAS No. : 3568-51-2  
 Formula: C<sub>6</sub>H<sub>10</sub>Cl<sub>3</sub>O<sub>4</sub>P  
 Molecular Weight: 283.47  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	AI 3-22016 is an insecticide.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5277 mL	17.6386 mL	35.2771 mL
5 mM	0.7055 mL	3.5277 mL	7.0554 mL
10 mM	0.3528 mL	1.7639 mL	3.5277 mL
50 mM	0.0706 mL	0.3528 mL	0.7055 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481