

N3-C5-NHS ester

Chemical Properties

CAS No. : 866363-70-4


Formula: C₁₀H₁₄N₄O₄

Molecular Weight: 254.24

Keep away from direct sunlight

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N3-C5-NHS ester, a noncleavable ADC linker, can be used in the synthesis of antibody-drug conjugates (ADCs).
Targets(IC50)	ADC Linker, PROTAC Linker

Solubility Information

Solubility	DMSO: 250 mg/mL (983.32 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4 mg/mL (15.73 mM), Sonication is recommended. 10% DMSO+90% Saline: 10 mg/mL (39.33 mM), Solution. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9333 mL	19.6665 mL	39.3329 mL
5 mM	0.7867 mL	3.9333 mL	7.8666 mL
10 mM	0.3933 mL	1.9666 mL	3.9333 mL
50 mM	0.0787 mL	0.3933 mL	0.7867 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Pyrrlobenzodiazepine Dimer Antibody-Drug Conjugates: Synthesis and Evaluation of Noncleavable Drug-Linkers. J Med Chem. 2017 Dec 14;60(23):9490-9507.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481