


pyridazin-4-ol

Chemical Properties

CAS No. :	17417-57-1
Formula:	C ₄ H ₄ N ₂ O
Molecular Weight:	96.08
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	pyridazin-4-ol, with CAS No. 17417-57-1, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. pyridazin-4-ol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	10.408 mL	52.040 mL	104.0799 mL
5 mM	2.0816 mL	10.408 mL	20.816 mL
10 mM	1.0408 mL	5.204 mL	10.408 mL
50 mM	0.2082 mL	1.0408 mL	2.0816 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481