


L 702958

## Chemical Properties

CAS No. : 136078-58-5  
 Formula: C<sub>25</sub>H<sub>28</sub>ClN<sub>3</sub>O<sub>4</sub>S  
 Molecular Weight: 502.03  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

|               |   |
|---------------|---|
| Description   | L 702958 is one kind of spirobenzopyran piperidine. |
| Targets(IC50) | Others  |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg       |
|-------|-----------|-----------|------------|
| 1 mM  | 1.9919 mL | 9.9596 mL | 19.9191 mL |
| 5 mM  | 0.3984 mL | 1.9919 mL | 3.9838 mL  |
| 10 mM | 0.1992 mL | 0.996 mL  | 1.9919 mL  |
| 50 mM | 0.0398 mL | 0.1992 mL | 0.3984 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Lynch JJ Jr, Wallace AA, Stupienski RF 3rd, Baskin EP, Beare CM, Appleby SD, Salata JJ, Jurkiewicz NK, Sanguinetti MC, Stein RB, et al. Cardiac electrophysiologic and antiarrhythmic actions of two long-acting spirobenzopyran piperidine class III agents, L-702,958 and L-706,000 [MK-499]. J Pharmacol Exp Ther. 1994 May;269(2):541-54. PubMed PMID: 7802864.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481