


3-Amino-2,6-piperidinedione

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 2353-44-8 |
| Formula: | C ₅ H ₈ N ₂ O ₂ |
| Molecular Weight: | 128.131 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|-------------|---|
| Description | 3-Amino-2,6-piperidinedione, with CAS No. 2353-44-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3-Amino-2,6-piperidinedione provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research. |
|-------------|---|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 7.8046 mL | 39.0229 mL | 78.0457 mL |
| 5 mM | 1.5609 mL | 7.8046 mL | 15.6091 mL |
| 10 mM | 0.7805 mL | 3.9023 mL | 7.8046 mL |
| 50 mM | 0.1561 mL | 0.7805 mL | 1.5609 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481