


Dichotomin

Chemical Properties

CAS No. :	53093-47-3
Formula:	C ₅₇ H ₉₄ O ₂₆
Molecular Weight:	1195.353
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dichotomin is a natural product of Paris, Liliaceae. The catalog number is TMA2055 and the CAS number is 53093-47-3. Dichotomin can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8366 mL	4.1829 mL	8.3658 mL
5 mM	0.1673 mL	0.8366 mL	1.6732 mL
10 mM	0.0837 mL	0.4183 mL	0.8366 mL
50 mM	0.0167 mL	0.0837 mL	0.1673 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481