

## 3'-Hydroxygynuramide II

## Chemical Properties


CAS No. : 1401093-57-9

Formula: C<sub>42</sub>H<sub>83</sub>N<sub>6</sub>O<sub>6</sub>

Molecular Weight: 698.11

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3'-Hydroxygynuramide II is a natural product for research related to life sciences. The catalog number is TN5821 and the CAS number is 1401093-57-9.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4324 mL	7.1622 mL	14.3244 mL
5 mM	0.2865 mL	1.4324 mL	2.8649 mL
10 mM	0.1432 mL	0.7162 mL	1.4324 mL
50 mM	0.0286 mL	0.1432 mL	0.2865 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481