


G-33378

## Chemical Properties

CAS No. : 1174-98-7  
 Formula: C<sub>23</sub>H<sub>20</sub>N<sub>2</sub>O<sub>3</sub>S  
 Molecular Weight: 404.48  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	G-33378 is a metabolite of sulfapyrazone that inhibits human platelet cyclo-oxygenase activity.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4723 mL	12.3616 mL	24.7231 mL
5 mM	0.4945 mL	2.4723 mL	4.9446 mL
10 mM	0.2472 mL	1.2362 mL	2.4723 mL
50 mM	0.0494 mL	0.2472 mL	0.4945 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481