


1-Methyl-3-phenyl-1,2,3,6-tetrahydropyridine

Chemical Properties

CAS No. :	119375-00-7
Formula:	C ₁₂ H ₁₅ N
Molecular Weight:	173.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Methyl-3-phenyl-1,2,3,6-tetrahydropyridine is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.772 mL	28.860 mL	57.7201 mL
5 mM	1.1544 mL	5.772 mL	11.544 mL
10 mM	0.5772 mL	2.886 mL	5.772 mL
50 mM	0.1154 mL	0.5772 mL	1.1544 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481