

pGlu-K(-NBD)-Q-R-L-G-N-Q-W-A-V-G-H-L-M-N

Chemical Properties

CAS No. :


Formula: C83H120N28O23S

Molecular Weight: 1910.08

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	pGlu-K(-NBD)-Q-R-L-G-N-Q-W-A-V-G-H-L-M-N is a peptides
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5235 mL	2.6177 mL	5.2354 mL
5 mM	0.1047 mL	0.5235 mL	1.0471 mL
10 mM	0.0524 mL	0.2618 mL	0.5235 mL
50 mM	0.0105 mL	0.0524 mL	0.1047 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481