

Ebracteolatanolide A

Chemical Properties


CAS No. :

Formula: C₂₀H₂₈O₄

Molecular Weight: 332.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | Ebracteolatanolide A is a useful organic compound for research related to life sciences and the catalog number is T125350. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 3.0081 mL | 15.0403 mL | 30.0806 mL |
| 5 mM | 0.6016 mL | 3.0081 mL | 6.0161 mL |
| 10 mM | 0.3008 mL | 1.504 mL | 3.0081 mL |
| 50 mM | 0.0602 mL | 0.3008 mL | 0.6016 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481