


TNFalpha-IN-S10

Chemical Properties

CAS No. :	920116-81-0
Formula:	C ₂₁ H ₁₄ N ₂ O ₃ S
Molecular Weight:	374.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	TNFalpha-IN-S10 is an inhibitor of TNF-α.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6709 mL	13.3543 mL	26.7087 mL
5 mM	0.5342 mL	2.6709 mL	5.3417 mL
10 mM	0.2671 mL	1.3354 mL	2.6709 mL
50 mM	0.0534 mL	0.2671 mL	0.5342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481