


Aminoacyl tRNA synthetase-IN-1

Chemical Properties

CAS No. :	219931-45-0
Formula:	C ₁₆ H ₂₅ N ₇ O ₇ S
Molecular Weight:	459.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Aminoacyl tRNA synthetase-IN-1 is an inhibitor of bacterial aminoacyl tRNA synthetase (aaRS).
Targets(IC50)	Others,Antibacterial

Solubility Information

Solubility	DMSO: 65 mg/mL (141.46 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 1.67 mg/mL (3.63 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1764 mL	10.8819 mL	21.7637 mL
5 mM	0.4353 mL	2.1764 mL	4.3527 mL
10 mM	0.2176 mL	1.0882 mL	2.1764 mL
50 mM	0.0435 mL	0.2176 mL	0.4353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Vondenhoff GH, et al. Microcin C and albomycin analogues with aryl-tetrazole substituents as nucleobase isosters are selective inhibitors of bacterial aminoacyl tRNA synthetases but lack efficient uptake. *Chembiochem*. 2012 Sep 3;13(13):1959-69.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481