

STAT3-IN-A69

## Chemical Properties


CAS No. : 1164546-70-6

Formula: C<sub>23</sub>H<sub>16</sub>N<sub>2</sub>O<sub>4</sub>

Molecular Weight: 384.38

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	STAT3-IN-A69 is an inhibitor targeting the DNA-binding domain of STAT3, suppressing tumor growth, metastasis and STAT3 target gene expression in vivo.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6016 mL	13.008 mL	26.0159 mL
5 mM	0.5203 mL	2.6016 mL	5.2032 mL
10 mM	0.2602 mL	1.3008 mL	2.6016 mL
50 mM	0.052 mL	0.2602 mL	0.5203 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481