

## 4-Cyano-4-phenylcyclohexanone

## Chemical Properties

CAS No. :	25115-74-6
Formula:	C13H13NO
Molecular Weight:	199.248
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	4-Cyano-4-phenylcyclohexanone, with CAS No. 25115-74-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Cyano-4-phenylcyclohexanone provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0188 mL	25.0941 mL	50.1882 mL
5 mM	1.0038 mL	5.0188 mL	10.0376 mL
10 mM	0.5019 mL	2.5094 mL	5.0188 mL
50 mM	0.1004 mL	0.5019 mL	1.0038 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481