


Pimprinine

Chemical Properties

CAS No. :	13640-26-1
Formula:	C ₁₂ H ₁₀ N ₂ O
Molecular Weight:	198.22
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pimprinine is an alkaloid. It is a potent inhibitor of MAO (monoamine oxidase). More recently, pimprinine showed promising anticonvulsant activity in electric seizure threshold.
Targets(IC50)	Others, Monoamine Oxidase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0449 mL	25.2245 mL	50.449 mL
5 mM	1.009 mL	5.0449 mL	10.0898 mL
10 mM	0.5045 mL	2.5224 mL	5.0449 mL
50 mM	0.1009 mL	0.5045 mL	1.009 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481