


Filic-3-en-25-al

Chemical Properties

CAS No. :	147850-78-0
Formula:	C ₃₀ H ₄₈ O
Molecular Weight:	424.713
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Filic-3-en-25-al is a natural product for research related to life sciences. The catalog number is TN4057 and the CAS number is 147850-78-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3545 mL	11.7727 mL	23.5455 mL
5 mM	0.4709 mL	2.3545 mL	4.7091 mL
10 mM	0.2355 mL	1.1773 mL	2.3545 mL
50 mM	0.0471 mL	0.2355 mL	0.4709 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481