


Compound X3805

Chemical Properties

CAS No. :	1426129-50-1
Formula:	C ₂₀ H ₂₅ NO ₃
Molecular Weight:	327.424
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound X3805, with CAS No. 1426129-50-1, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Compound X3805 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0542 mL	15.2709 mL	30.5418 mL
5 mM	0.6108 mL	3.0542 mL	6.1084 mL
10 mM	0.3054 mL	1.5271 mL	3.0542 mL
50 mM	0.0611 mL	0.3054 mL	0.6108 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481