

HAEGT

Chemical Properties

CAS No. : 852155-81-8


Formula: C₂₀H₃₁N₇O₉

Molecular Weight: 513.5

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	HAEGT is the N-terminal 1-5 residues of the GLP-1 peptide (glucagon-like peptide-1).
Targets(IC ₅₀)	Proteasome,Glucagon Receptor

Solubility Information

Solubility	H ₂ O: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9474 mL	9.7371 mL	19.4742 mL
5 mM	0.3895 mL	1.9474 mL	3.8948 mL
10 mM	0.1947 mL	0.9737 mL	1.9474 mL
50 mM	0.0389 mL	0.1947 mL	0.3895 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481