


## 2-Benzylthioadenosine

## Chemical Properties

CAS No. :	43157-48-8
Formula:	C17H19N5O4S
Molecular Weight:	389.43
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Benzylthioadenosine is a Nucleoside Derivative - 2-Modified purine nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5679 mL	12.8393 mL	25.6786 mL
5 mM	0.5136 mL	2.5679 mL	5.1357 mL
10 mM	0.2568 mL	1.2839 mL	2.5679 mL
50 mM	0.0514 mL	0.2568 mL	0.5136 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481