


Butyric acid, 2-ethyl-4-hydroxy-, (5-nitrofurfurylidene)hydrazide

Chemical Properties

CAS No. :	91338-37-3
Formula:	C ₁₁ H ₁₅ N ₃ O ₅
Molecular Weight:	269.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Butyric acid, 2-ethyl-4-hydroxy-, (5-nitrofurfurylidene)hydrazide is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.714 mL	18.5701 mL	37.1402 mL
5 mM	0.7428 mL	3.714 mL	7.428 mL
10 mM	0.3714 mL	1.857 mL	3.714 mL
50 mM	0.0743 mL	0.3714 mL	0.7428 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481