

Compound N039-0044

## Chemical Properties


CAS No. :

Formula: C37H60O13

Molecular Weight: 712.874

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N039-0044 is a useful organic compound for research related to life sciences and the catalog number is T131266.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4028 mL	7.0139 mL	14.0278 mL
5 mM	0.2806 mL	1.4028 mL	2.8056 mL
10 mM	0.1403 mL	0.7014 mL	1.4028 mL
50 mM	0.0281 mL	0.1403 mL	0.2806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481