

Compound N045-0004

Chemical Properties


CAS No. :

Formula: C40H60N2O4

Molecular Weight: 632.93

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N045-0004 is a useful organic compound for research related to life sciences and the catalog number is T131255.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.580 mL	7.8998 mL	15.7995 mL
5 mM	0.316 mL	1.580 mL	3.1599 mL
10 mM	0.158 mL	0.790 mL	1.580 mL
50 mM	0.0316 mL	0.158 mL	0.316 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481