

Compound 0449-0145

## Chemical Properties


CAS No. :

Formula: C49H76O20

Molecular Weight: 985.127

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0449-0145 is a useful organic compound for research related to life sciences and the catalog number is T131643.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0151 mL	5.0755 mL	10.1509 mL
5 mM	0.203 mL	1.0151 mL	2.0302 mL
10 mM	0.1015 mL	0.5075 mL	1.0151 mL
50 mM	0.0203 mL	0.1015 mL	0.203 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481