


S-3304

## Chemical Properties

CAS No. : 193809-34-6  
 Formula: C<sub>24</sub>H<sub>20</sub>N<sub>2</sub>O<sub>4</sub>S<sub>2</sub>  
 Molecular Weight: 464.56  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	S-3304 is a n orally-agent agent with potential antineoplastic activity. S-3304 inhibits matrix metalloproteinases (MMPs), thereby inducing extracellular matrix degradation, and inhibiting angiogenesis, tumor growth and invasion, and metastasis.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1526 mL	10.7629 mL	21.5257 mL
5 mM	0.4305 mL	2.1526 mL	4.3051 mL
10 mM	0.2153 mL	1.0763 mL	2.1526 mL
50 mM	0.0431 mL	0.2153 mL	0.4305 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481