


## [Hyp3]-Bradykinin

## Chemical Properties

CAS No. :	37642-65-2
Formula:	C50H73N15O12
Molecular Weight:	1076.227
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Bradykinin antagonist.
Targets(IC50)	Bradykinin Receptor

## Solubility Information

Solubility	DMSO: Soluble ( $< 1$ mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9292 mL	4.6458 mL	9.2917 mL
5 mM	0.1858 mL	0.9292 mL	1.8583 mL
10 mM	0.0929 mL	0.4646 mL	0.9292 mL
50 mM	0.0186 mL	0.0929 mL	0.1858 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481