

Compound T128140

Chemical Properties


CAS No. :

Formula: C₂₂H₂₄O₁₃

Molecular Weight: 496.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

4-hydroxy-3-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-((4-hydroxy-3,5-dimethoxybenzoyloxy)methyl)oxan-2-yl]oxy)benzoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0144 mL	10.0721 mL	20.1442 mL
5 mM	0.4029 mL	2.0144 mL	4.0288 mL
10 mM	0.2014 mL	1.0072 mL	2.0144 mL
50 mM	0.0403 mL	0.2014 mL	0.4029 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481