

Tuberostemonine D

Chemical Properties

CAS No. :	1627827-84-2
Formula:	C ₂₂ H ₃₃ N ₄ O ₄
Molecular Weight:	375.50
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Tuberostemonine D is a useful organic compound for research related to life sciences and the catalog number is T125659.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6631 mL	13.3156 mL	26.6312 mL
5 mM	0.5326 mL	2.6631 mL	5.3262 mL
10 mM	0.2663 mL	1.3316 mL	2.6631 mL
50 mM	0.0533 mL	0.2663 mL	0.5326 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481