


ERR α antagonist-1

Chemical Properties

CAS No. :	1072145-33-5
Formula:	C ₂₁ H ₁₉ N ₃ S ₂
Molecular Weight:	377.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	ERR α antagonist-1 (ERR+/- antagonist-1) is a high-affinity, selective antagonist of the estrogen-related receptor α (ERR α). It effectively prevents the interaction of ERR α with both Proliferator-activated Receptor γ Coactivator-1 α (PGC-1 α) and PGC-1 β , displaying IC ₅₀ values of 170 nM and 180 nM, respectively. Notably, ERR α antagonist-1 does not interfere with the interactions involving ERR β or ERR γ and the PGC-1 α and PGC-1 β coactivators. Furthermore, it does not affect the interaction between either ER α or ER β and PGC-1 α or SRC-1.
Targets(IC ₅₀)	Estrogen Receptor/ERR, Estrogen/progestogen Receptor

Solubility Information

Solubility	DMSO: 8 mg/mL (21.19 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6488 mL	13.244 mL	26.488 mL
5 mM	0.5298 mL	2.6488 mL	5.2976 mL
10 mM	0.2649 mL	1.3244 mL	2.6488 mL
50 mM	0.053 mL	0.2649 mL	0.5298 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Michael J. Chisamore, et al. Identification of Small Molecule Estrogen-Related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481