

Ethyl ganoderate J

Chemical Properties


CAS No. :

Formula: C32H46O7

Molecular Weight: 542.713

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethyl ganoderate J is a useful organic compound for research related to life sciences and the catalog number is T126258.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8426 mL	9.213 mL	18.426 mL
5 mM	0.3685 mL	1.8426 mL	3.6852 mL
10 mM	0.1843 mL	0.9213 mL	1.8426 mL
50 mM	0.0369 mL	0.1843 mL	0.3685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481