

Deldeprevir Na

Chemical Properties

CAS No. :	1298053-61-8
Formula:	C ₄₅ H ₅₅ F ₂ N ₆ NaO ₈ S ₂
Molecular Weight:	933.08
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Deldeprevir Na is an Antiviral Polyprotein cleavage inhibitor. HCV NS3/4A inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0717 mL	5.3586 mL	10.7172 mL
5 mM	0.2143 mL	1.0717 mL	2.1434 mL
10 mM	0.1072 mL	0.5359 mL	1.0717 mL
50 mM	0.0214 mL	0.1072 mL	0.2143 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481