


SGC-AAK1-1N

Chemical Properties

CAS No. :	2242913-80-8
Formula:	C ₂₀ H ₂₂ N ₄ O ₃ S
Molecular Weight:	398.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	SGC-AAK1-1N is a potent and selective inhibitor of AAK1 (AP2 associated kinase 1), with an IC ₅₀ of 1.8 μM [1].
Targets(IC ₅₀)	Others,AAK1 (AP2 associated kinase 1)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5095 mL	12.5477 mL	25.0954 mL
5 mM	0.5019 mL	2.5095 mL	5.0191 mL
10 mM	0.251 mL	1.2548 mL	2.5095 mL
50 mM	0.0502 mL	0.251 mL	0.5019 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481