

PDE4B-IN-3

Chemical Properties

CAS No. :

Formula: C30H35N3O4S2

Molecular Weight: 565.75

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PDE4B-IN-3 is a potent PDE4B inhibitor (IC50: 0.94 μ M) with anti-inflammatory effects.
Targets(IC50)	Others,PDE

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7676 mL	8.8378 mL	17.6757 mL
5 mM	0.3535 mL	1.7676 mL	3.5351 mL
10 mM	0.1768 mL	0.8838 mL	1.7676 mL
50 mM	0.0354 mL	0.1768 mL	0.3535 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481