


1-Benzofuran-5-ylmethanol

Chemical Properties

CAS No. :	31823-05-9
Formula:	C ₉ H ₈ O ₂
Molecular Weight:	148.16
Storage:	Store at low temperature
	Pure form: -20°C for 3 years In solvent: -80°C for 1 year
	<small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	1-Benzofuran-5-ylmethanol, with CAS No. 31823-05-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1-Benzofuran-5-ylmethanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7495 mL	33.7473 mL	67.4946 mL
5 mM	1.3499 mL	6.7495 mL	13.4989 mL
10 mM	0.6749 mL	3.3747 mL	6.7495 mL
50 mM	0.135 mL	0.6749 mL	1.3499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481