


DS1

Chemical Properties

CAS No. : 372497-52-4
 Formula: C₁₈H₁₀Br₂ClN₃OS
 Molecular Weight: 511.62
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	DS1 is a GABAA receptor agonist which can stimulate gonadotropin subunit gene expression in mouse.
Targets(IC50)	Others,GABA Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9546 mL	9.7729 mL	19.5458 mL
5 mM	0.3909 mL	1.9546 mL	3.9092 mL
10 mM	0.1955 mL	0.9773 mL	1.9546 mL
50 mM	0.0391 mL	0.1955 mL	0.3909 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481