


Arenarine D

Chemical Properties

CAS No. :	123520-96-7
Formula:	C ₁₃ H ₁₀ N ₂ O ₂
Molecular Weight:	226.235
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Arenarine D is a useful organic compound for research related to life sciences. The catalog number is T126174 and the CAS number is 123520-96-7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4201 mL	22.1004 mL	44.2008 mL
5 mM	0.884 mL	4.4201 mL	8.8402 mL
10 mM	0.442 mL	2.210 mL	4.4201 mL
50 mM	0.0884 mL	0.442 mL	0.884 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481