


Acalisib

Chemical Properties

CAS No. :	870281-34-8
Formula:	C ₂₁ H ₁₆ N ₇ O
Molecular Weight:	401.4
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acalisib (CAL-120) (GS-9820) is a potent and selective inhibitor of PI3Kδ.
Targets(IC50)	PI3K

Solubility Information

Solubility	DMSO: 73 mg/mL (181.86 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (4.98 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4913 mL	12.4564 mL	24.9128 mL
5 mM	0.4983 mL	2.4913 mL	4.9826 mL
10 mM	0.2491 mL	1.2456 mL	2.4913 mL
50 mM	0.0498 mL	0.2491 mL	0.4983 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Shugg RP, et al. J Biol Chem. 2013 Dec 6;288(49):35346-57.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481