

Otenaproxesul

Chemical Properties


CAS No. : 1226895-20-0

Formula: C₂₁H₁₉N₃O₃S

Molecular Weight: 365.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Otenaproxesul (ATB 346) is a novel hydrogen sulphide-releasing derivative of naproxen with markedly reduced toxicity.
Targets(IC50)	Apoptosis, COX

Solubility Information

Solubility	DMSO: 50 mg/mL (136.82 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.47 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7364 mL	13.6818 mL	27.3635 mL
5 mM	0.5473 mL	2.7364 mL	5.4727 mL
10 mM	0.2736 mL	1.3682 mL	2.7364 mL
50 mM	0.0547 mL	0.2736 mL	0.5473 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wallace JL, et al. Br J Pharmacol. 2010 Mar;159(6):1236-46.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481