

Compound N005-0020

## Chemical Properties


CAS No. :

Formula: C<sub>35</sub>H<sub>32</sub>O<sub>15</sub>

Molecular Weight: 692.626

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N005-0020 is a useful organic compound for research related to life sciences and the catalog number is T131375.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4438 mL	7.2189 mL	14.4377 mL
5 mM	0.2888 mL	1.4438 mL	2.8875 mL
10 mM	0.1444 mL	0.7219 mL	1.4438 mL
50 mM	0.0289 mL	0.1444 mL	0.2888 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481