


Ginnol

Chemical Properties

CAS No. :	2606-50-0
Formula:	C ₂₉ H ₆₀ O
Molecular Weight:	424.79
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ginno is a natural product Ginkgo biloba L.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3541 mL	11.7705 mL	23.541 mL
5 mM	0.4708 mL	2.3541 mL	4.7082 mL
10 mM	0.2354 mL	1.1771 mL	2.3541 mL
50 mM	0.0471 mL	0.2354 mL	0.4708 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481