


## (3-Pyrimidin-5-ylphenyl)methanol

## Chemical Properties

CAS No. :	852180-75-7
Formula:	C <sub>11</sub> H <sub>10</sub> N <sub>2</sub> O
Molecular Weight:	186.21
Storage:	Store under nitrogen Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	(3-Pyrimidin-5-ylphenyl)methanol, with CAS No. 852180-75-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (3-Pyrimidin-5-ylphenyl)methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3703 mL	26.8514 mL	53.7028 mL
5 mM	1.0741 mL	5.3703 mL	10.7406 mL
10 mM	0.537 mL	2.6851 mL	5.3703 mL
50 mM	0.1074 mL	0.537 mL	1.0741 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481