

## 3,8a-dihydroxy-3,8-dimethyl-5-(propan-2-ylidene)-1,2,3,3a,4,5,6,8a-octahydroazulen-6-one

## Chemical Properties


CAS No. :

Formula: C<sub>15</sub>H<sub>22</sub>O<sub>3</sub>

Molecular Weight: 250.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3,8a-dihydroxy-3,8-dimethyl-5-(propan-2-ylidene)-1,2,3,3a,4,5,6,8a-octahydroazulen-6-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9946 mL	19.9728 mL	39.9457 mL
5 mM	0.7989 mL	3.9946 mL	7.9891 mL
10 mM	0.3995 mL	1.9973 mL	3.9946 mL
50 mM	0.0799 mL	0.3995 mL	0.7989 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481