

Pd-C-II

Chemical Properties

CAS No. :

Formula: C₁₉H₂₀O₆

Molecular Weight: 344.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pd-C-II has anti-inflammation activity, it can inhibit TNF- α production and iNOS protein expression and inhibit COX-2 protein expression in LPS-stimulated RAW 264.7 cells. Pd-C-II inhibits anaphylactic mediator release from purified mast cells induced by concanavalin A with phosphatidylserine, the IC ₅₀ value of 79 microM.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9036 mL	14.518 mL	29.036 mL
5 mM	0.5807 mL	2.9036 mL	5.8072 mL
10 mM	0.2904 mL	1.4518 mL	2.9036 mL
50 mM	0.0581 mL	0.2904 mL	0.5807 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481