

Iodo-PEG4-N3

Chemical Properties

CAS No. : 936917-36-1


Formula: C₈H₁₆IN₃O₃

Molecular Weight: 329.14

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Iodo-PEG4-N3 is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, comprising two ligands joined by a linker—one binding an E3 ubiquitin ligase and the other the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0382 mL	15.1911 mL	30.3822 mL
5 mM	0.6076 mL	3.0382 mL	6.0764 mL
10 mM	0.3038 mL	1.5191 mL	3.0382 mL
50 mM	0.0608 mL	0.3038 mL	0.6076 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481