


3,7,16-Trihydroxystigmast-5-ene

Chemical Properties

CAS No. :	289056-24-2
Formula:	C ₂₉ H ₅₀ O ₃
Molecular Weight:	446.716
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,7,16-Trihydroxystigmast-5-ene is a natural product for research related to life sciences. The catalog number is TN2908 and the CAS number is 289056-24-2.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2385 mL	11.1927 mL	22.3854 mL
5 mM	0.4477 mL	2.2385 mL	4.4771 mL
10 mM	0.2239 mL	1.1193 mL	2.2385 mL
50 mM	0.0448 mL	0.2239 mL	0.4477 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481