


## 2-PHENOXYBENZONITRILE

## Chemical Properties

CAS No. :	6476-32-0
Formula:	C <sub>10</sub> H <sub>9</sub> NO
Molecular Weight:	159.18
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-PHENOXYBENZONITRILE, with CAS No. 6476-32-0, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-PHENOXYBENZONITRILE provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.2822 mL	31.411 mL	62.822 mL
5 mM	1.2564 mL	6.2822 mL	12.5644 mL
10 mM	0.6282 mL	3.1411 mL	6.2822 mL
50 mM	0.1256 mL	0.6282 mL	1.2564 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481