

Compound T129036

Chemical Properties


CAS No. :

Formula: C53H84O23

Molecular Weight: 1089.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3S,4S,5R,6R)-6-(((3S,6aR,6bS,8aS,12aS,14bR)-4,4,6a,6b,11,11,14b-heptamethyl-8a-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)carbonyl)-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydricen-3-yl]oxy)-3-hydroxy-5-(((2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)-4-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)oxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9181 mL	4.5904 mL	9.1808 mL
5 mM	0.1836 mL	0.9181 mL	1.8362 mL
10 mM	0.0918 mL	0.459 mL	0.9181 mL
50 mM	0.0184 mL	0.0918 mL	0.1836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481