


Mexrenoate potassium

Chemical Properties

CAS No. :	41020-67-1
Formula:	C ₂₄ H ₃₃ KO ₆
Molecular Weight:	456.61
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mexrenoate potassium is a steroidal aldosterone antagonist and antihypertensive.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1901 mL	10.9503 mL	21.9005 mL
5 mM	0.438 mL	2.1901 mL	4.3801 mL
10 mM	0.219 mL	1.095 mL	2.1901 mL
50 mM	0.0438 mL	0.219 mL	0.438 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481