

Compound T130403

Chemical Properties

CAS No. :

Formula: C₃₂H₅₆O₁₅

Molecular Weight: 680.79

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3-methyl-1-[[3-methyl-1-({3-methyl-1-oxo-1-[(2,3,4,5,6-pentahydroxyhexyl)oxy]pentan-2-yl]oxy)-1-oxopentan-2-yl]oxy]-1-oxopentan-2-yl 2-(acetyloxy)-3-methylpentanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4689 mL	7.3444 mL	14.6888 mL
5 mM	0.2938 mL	1.4689 mL	2.9378 mL
10 mM	0.1469 mL	0.7344 mL	1.4689 mL
50 mM	0.0294 mL	0.1469 mL	0.2938 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481